


S(+)-MDO-NPA HCl

Chemical Properties

CAS No. :	113678-73-2
Formula:	C ₂₀ H ₂₂ ClNO ₂
Molecular Weight:	343.85
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	S(+)-MDO-NPA HCl is an inhibitor of dopamine with oral activity, selective for rat limbic system.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9082 mL	14.5412 mL	29.0824 mL
5 mM	0.5816 mL	2.9082 mL	5.8165 mL
10 mM	0.2908 mL	1.4541 mL	2.9082 mL
50 mM	0.0582 mL	0.2908 mL	0.5816 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481