

## TCO-PEG3-oxyamine

## Chemical Properties

CAS No. :


Formula: C19H35N3O7

Molecular Weight: 417.5

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	TCO-PEG3-oxyamine is a PEG-based linker for PROTACs that joins two essential ligands, facilitating selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two distinct ligands linked together: one binds to an E3 ubiquitin ligase and the other targets a specific protein. They leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3952 mL	11.976 mL	23.9521 mL
5 mM	0.479 mL	2.3952 mL	4.7904 mL
10 mM	0.2395 mL	1.1976 mL	2.3952 mL
50 mM	0.0479 mL	0.2395 mL	0.479 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481