


Ono 11119

Chemical Properties

CAS No. : 82402-68-4
 Formula: C₂₆H₄₇NO₃
 Molecular Weight: 421.66
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ono 11119 is a thromboxane A2 antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3716 mL	11.8579 mL	23.7158 mL
5 mM	0.4743 mL	2.3716 mL	4.7432 mL
10 mM	0.2372 mL	1.1858 mL	2.3716 mL
50 mM	0.0474 mL	0.2372 mL	0.4743 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481