

## Ethyl 1-benzofuran-3-ylacetate

## Chemical Properties


CAS No. : 82156-58-9

Formula: C<sub>12</sub>H<sub>12</sub>O<sub>3</sub>

Molecular Weight: 204.22

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

## Description

Ethyl 1-benzofuran-3-ylacetate, with CAS No. 82156-58-9, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. Ethyl 1-benzofuran-3-ylacetate provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8967 mL	24.4834 mL	48.9668 mL
5 mM	0.9793 mL	4.8967 mL	9.7934 mL
10 mM	0.4897 mL	2.4483 mL	4.8967 mL
50 mM	0.0979 mL	0.4897 mL	0.9793 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481