


Ulonivirine

Chemical Properties

CAS No. :	1591823-76-5
Formula:	C ₁₈ H ₈ ClF ₆ N ₅ O ₃
Molecular Weight:	491.73
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ulonivirine is an orally active non-nucleoside reverse transcriptase inhibitor that exhibits high antiviral activity and can be used to study HIV-1 infection.
Targets(IC50)	Others,HIV Protease,Reverse Transcriptase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0336 mL	10.1682 mL	20.3364 mL
5 mM	0.4067 mL	2.0336 mL	4.0673 mL
10 mM	0.2034 mL	1.0168 mL	2.0336 mL
50 mM	0.0407 mL	0.2034 mL	0.4067 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481