

m-PEG-NH2 (MW 1000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	m-PEG-NH2 (MW 1000) is a PEG-based PROTAC linker. m-PEG-amine (MW 1000) can be used to synthesize PROTAC.
Targets(IC50)	PROTAC Linker

Solubility Information

Solubility	H2O: 200 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481