


2',4'-Dihydroxy-3,7':4,8'-diepoxylign-7-ene

Chemical Properties

CAS No. : 666250-52-8
 Formula: C₁₈H₁₈O₄
 Molecular Weight: 298.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	2',4'-Dihydroxy-3,7':4,8'-diepoxylign-7-ene is a natural compound used in life sciences research, catalog number TN2716, CAS number [666250-52-8].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.352 mL	16.760 mL	33.5199 mL
5 mM	0.6704 mL	3.352 mL	6.704 mL
10 mM	0.3352 mL	1.676 mL	3.352 mL
50 mM	0.067 mL	0.3352 mL	0.6704 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481