

Azido-PEG4-azide

Chemical Properties

CAS No. : 182760-73-2


Formula: C10H20N6O4

Molecular Weight: 288.3

Keep away from direct sunlight

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azido-PEG4-azide is a PEG-based PROTAC linker used for synthesizing PROTACs.
Targets(IC50)	PROTAC Linker
In vitro	PROTACs comprise two distinct ligands joined by a linker: one ligand targets an E3 ubiquitin ligase, while the other binds to the target protein [1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4686 mL	17.343 mL	34.6861 mL
5 mM	0.6937 mL	3.4686 mL	6.9372 mL
10 mM	0.3469 mL	1.7343 mL	3.4686 mL
50 mM	0.0694 mL	0.3469 mL	0.6937 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481