


## Dehydrodynerigenin digitaloside

## Chemical Properties

CAS No. :	52628-62-3
Formula:	C30H42O8
Molecular Weight:	530.65
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dehydrodynerigenin digitaloside is a natural product for research related to life sciences. The catalog number is TN3791 and the CAS number is 52628-62-3.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8845 mL	9.4224 mL	18.8448 mL
5 mM	0.3769 mL	1.8845 mL	3.769 mL
10 mM	0.1884 mL	0.9422 mL	1.8845 mL
50 mM	0.0377 mL	0.1884 mL	0.3769 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481