


## CASP3 Activator 1541B

## Chemical Properties

CAS No. :	1170689-68-5
Formula:	C <sub>23</sub> H <sub>15</sub> N <sub>3</sub> O <sub>4</sub>
Molecular Weight:	397.38
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	CASP3 Activator 1541B is an activator of the proenzyme forms of caspase-3, and is an inducer of cell death.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5165 mL	12.5824 mL	25.1648 mL
5 mM	0.5033 mL	2.5165 mL	5.033 mL
10 mM	0.2516 mL	1.2582 mL	2.5165 mL
50 mM	0.0503 mL	0.2516 mL	0.5033 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481