


Proglumide sodium

Chemical Properties

CAS No. :	99247-33-3
Formula:	C ₁₈ H ₂₆ N ₂ NaO ₄
Molecular Weight:	357.4
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Non-selective cholecystokinin (CCK) antagonist
Targets(IC50)	Others,Cholecystokinin Receptor

Solubility Information

Solubility	H ₂ O: <35.64 mg/mL, Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.798 mL	13.9899 mL	27.9799 mL
5 mM	0.5596 mL	2.798 mL	5.596 mL
10 mM	0.2798 mL	1.399 mL	2.798 mL
50 mM	0.056 mL	0.2798 mL	0.5596 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481