

C.I. 59830

## Chemical Properties

CAS No. : 25704-81-8


Formula: C<sub>36</sub>H<sub>18</sub>Br<sub>2</sub>O<sub>4</sub>

Molecular Weight: 674.344

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	C.I. 59830 is a green dye.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4829 mL	7.4147 mL	14.8293 mL
5 mM	0.2966 mL	1.4829 mL	2.9659 mL
10 mM	0.1483 mL	0.7415 mL	1.4829 mL
50 mM	0.0297 mL	0.1483 mL	0.2966 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481