


m-Anisaldehyde

Chemical Properties

CAS No. :	591-31-1
Formula:	C ₈ H ₈ O ₂
Molecular Weight:	136.15
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	m-Anisaldehyde exhibits binding activity toward aldo-keto reductase family 1 member A1 and is widely used in biochemical experiments and drug synthesis research.
Targets(IC50)	Others,Endogenous Metabolite

Solubility Information

Solubility	DMSO: 80 mg/mL (587.59 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 3.3 mg/mL (24.24 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.3448 mL	36.7242 mL	73.4484 mL
5 mM	1.469 mL	7.3448 mL	14.6897 mL
10 mM	0.7345 mL	3.6724 mL	7.3448 mL
50 mM	0.1469 mL	0.7345 mL	1.469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Alan CuthbertsonMagne Solbakken. Imaging and radiotherapy methods for tumour stem cells. WO2013048811 A1

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481