

AM-1430

## Chemical Properties


CAS No. :

Formula: C<sub>23</sub>H<sub>19</sub>N<sub>8</sub>O

Molecular Weight: 442.46

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	AM-1430 is a potent and selective PI3K $\delta$ inhibitor. AM-1430 has cellular potency: pAKT IC <sub>50</sub> = 0.8 nM; HWB,u (pAKT) IC <sub>50</sub> = 7.7 nM. Rat KLH: IgG ED <sub>50</sub> = 0.016 mg/kg; IhM ED <sub>50</sub> = 0.015 mg/kg. Rat PK Clu = 1.7 L/hr/kg.
Targets(IC <sub>50</sub> )	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2601 mL	11.3005 mL	22.6009 mL
5 mM	0.452 mL	2.2601 mL	4.5202 mL
10 mM	0.226 mL	1.130 mL	2.2601 mL
50 mM	0.0452 mL	0.226 mL	0.452 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Gonzalez-Lopez de Turiso F, Hao X, Shin Y, et al. Discovery and in Vivo Evaluation of the Potent and Selective PI3K $\delta$  Inhibitors 2-((1S)-1-((6-Amino-5-cyano-4-pyrimidinyl)amino)ethyl)-6-fluoro-N-methyl-3-(2-pyridinyl)-4-quinolinecarboxamide (AM-0687) and 2-((1S)-1-((6-Amino-5-cyano-4-pyrimidinyl)amino)ethyl)-5-fluoro-N-methyl-3-(2-pyridinyl)-4-quinolinecarboxamide (AM-1430). J Med Chem. 2016 Jul 14. [Epub ahead of print] PubMed PMID: 27411843.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481