

3-(4-Ethylbenzoyl)propionic acid

Chemical Properties


CAS No. : 49594-75-4

Formula: C₁₂H₁₄O₃

Molecular Weight: 206.24

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-(4-Ethylbenzoyl)propionic acid, with CAS No. 49594-75-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 3-(4-Ethylbenzoyl)propionic acid provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.8487 mL	24.2436 mL	48.4872 mL
5 mM	0.9697 mL	4.8487 mL	9.6974 mL
10 mM	0.4849 mL	2.4244 mL	4.8487 mL
50 mM	0.097 mL	0.4849 mL	0.9697 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481