


Ethyl 3,5-di-O-caffeoylquininate

Chemical Properties

CAS No. : 143051-74-5
 Formula: C₂₇H₂₈O₁₂
 Molecular Weight: 544.509
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ethyl 3,5-di-O-caffeoylquininate is a natural product for research related to life sciences. The catalog number is TN5844 and the CAS number is 143051-74-5.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8365 mL	9.1826 mL	18.3651 mL
5 mM	0.3673 mL	1.8365 mL	3.673 mL
10 mM	0.1837 mL	0.9183 mL	1.8365 mL
50 mM	0.0367 mL	0.1837 mL	0.3673 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481