

Lig3

Chemical Properties


CAS No. : 1448034-06-7

Formula: C₂₀H₁₇ClN₄O₃

Molecular Weight: 396.83

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lig3 is an inhibitor against high-risk strains of human papillomavirus (HPV) E6 oncoprotein.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.520 mL	12.5999 mL	25.1997 mL
5 mM	0.504 mL	2.520 mL	5.0399 mL
10 mM	0.252 mL	1.260 mL	2.520 mL
50 mM	0.0504 mL	0.252 mL	0.504 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481