


## Cannabicitran

## Chemical Properties

CAS No. :	31508-71-1
Formula:	C <sub>21</sub> H <sub>30</sub> O <sub>2</sub>
Molecular Weight:	314.46
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Cannabicitran is a cannabinoid that can decrease rabbits intraocular pressure [1] [2].
Targets(IC50)	Cannabinoid Receptor,Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1801 mL	15.9003 mL	31.8005 mL
5 mM	0.636 mL	3.1801 mL	6.3601 mL
10 mM	0.318 mL	1.590 mL	3.1801 mL
50 mM	0.0636 mL	0.318 mL	0.636 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481