

2,2-Dimethyl-2H-chromene-6-carbonitrile

Chemical Properties

CAS No. :	33143-29-2
Formula:	C ₁₂ H ₁₁ NO
Molecular Weight:	185.222
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

Biological Description

Description	2,2-Dimethyl-2H-chromene-6-carbonitrile, with CAS No. 33143-29-2, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. 2,2-Dimethyl-2H-chromene-6-carbonitrile provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.399 mL	26.9949 mL	53.9898 mL
5 mM	1.0798 mL	5.399 mL	10.798 mL
10 mM	0.5399 mL	2.6995 mL	5.399 mL
50 mM	0.108 mL	0.5399 mL	1.0798 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481