


Dihydrosphingosine

Chemical Properties

CAS No. :	3102-56-5
Formula:	C ₁₈ H ₃₉ NO ₂
Molecular Weight:	301.51
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Dihydrosphingosine is a Protein kinase C inhibitor.
Targets(IC50)	Others,Phospholipase,PKC

Solubility Information

Solubility	DMSO: 100 mM,Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3166 mL	16.5832 mL	33.1664 mL
5 mM	0.6633 mL	3.3166 mL	6.6333 mL
10 mM	0.3317 mL	1.6583 mL	3.3166 mL
50 mM	0.0663 mL	0.3317 mL	0.6633 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481