


## (3,5-Dimethyl-4-isoxazolyl)methanol

## Chemical Properties

CAS No. :	19788-36-4
Formula:	C <sub>6</sub> H <sub>9</sub> NO <sub>2</sub>
Molecular Weight:	127.14
Storage:	Store at low temperature Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	(3,5-Dimethyl-4-isoxazolyl)methanol, with CAS No. 19788-36-4, is a fragment molecule that serves as an important scaffold for molecular linking, expansion, and modification. (3,5-Dimethyl-4-isoxazolyl)methanol provides a structural basis and research tool for the design and screening of novel drug candidates, and is commonly used in drug discovery, drug synthesis, and related research.
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.8653 mL	39.3267 mL	78.6535 mL
5 mM	1.5731 mL	7.8653 mL	15.7307 mL
10 mM	0.7865 mL	3.9327 mL	7.8653 mL
50 mM	0.1573 mL	0.7865 mL	1.5731 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481