

c-Kit-IN-3 maleate

Chemical Properties


CAS No. :

Formula: C₃₀H₂₄ClF₃N₂O₈

Molecular Weight: 632.97

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	c-Kit-IN-3 maleate is a potent and selective c-KIT kinase inhibitor (IC ₅₀ s: 4 nM, 8 nM for c-KIT wt, and c-KIT T670I).
Targets(IC ₅₀)	Others
In vitro	c-Kit-IN-3 maleate (Compound 18; 0.1-10 μM; 6 days; primary GIST patient cells) exhibits dose-dependent antiproliferative effects. c-Kit-IN-3 (0.01-1 μM; 24 hours; GIST-T1, GIST-T1-T670I, and GIST-5R cells) treatment induces dose-dependent cell apoptotic death. c-Kit-IN-3 (0.01-1 μM; 24 hours; GIST-T1, GIST-T1-T670I, and GIST-5R cells) treatment arrests the cell cycle into the G ₀ /G ₁ phase in all of these three cell lines. c-Kit-IN-3 (0-1 μM; 24 hours) blocks the autophosphorylation of c-KIT pY703, pY719, and pY823 in GIST-T1, GIST-T1-T670I, and GIST-5R, respectively, cells at a concentration of 30 nM and inhibits the downstream signaling mediators pAKT (T308/S473), pS6 (S235/236), and pERK (T202/204). c-Kit-IN-3 potently inhibits the activity of CSF1R (IC ₅₀ : 18 nM), PDGFRα (IC ₅₀ : 25 nM), RET (IC ₅₀ : 34 nM), and it relatively less potently inhibits DDR1 (IC ₅₀ : 135 nM), FLT4 (IC ₅₀ : 121 nM), and PDGFRβ (IC ₅₀ : 97 nM). c-Kit-IN-3 (0.006 μM-1.37 μM) potently inhibits the growth of c-KIT-dependent GIST cancer cells, such as GIST-T1 (IC ₅₀ : 0.006 μM); GIST-882 (IC ₅₀ : 0.013 μM); GIST-T1-T670I (IC ₅₀ : 0.011 μM); GIST-5R (IC ₅₀ : 0.073 μM); GIST-48B (IC ₅₀ : 1.37 μM), respectively.
In vivo	c-Kit-IN-3 maleate (oral gavage; 20-100 mg/kg/day; for 11 days; female BALB/C-nu mice) treatment dose-dependently inhibits the tumor progression. c-Kit-IN-3 (1 mg/kg iv for mice, rats, dog; 10 mg/kg p.o. for mice, rats; and 5 mg/kg p.o. for dog) has T _{1/2} of 4.5 h, 6.4 h, 19.4 h for mice, rats and dogs, respectively. And it possesses acceptable bioavailability in mice (F = 43%), rats (F = 50%), and dogs (F = 81%).

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5799 mL	7.8993 mL	15.7985 mL
5 mM	0.316 mL	1.5799 mL	3.1597 mL
10 mM	0.158 mL	0.7899 mL	1.5799 mL
50 mM	0.0316 mL	0.158 mL	0.316 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wu Y, et al. Discovery of 2-(4-Chloro-3-(trifluoromethyl)phenyl)-N-(4-((6,7-dimethoxyquinolin-4-yl)oxy)phenyl)acetamide (CHMFL-KIT-64) as a Novel Orally Available Potent Inhibitor against Broad-Spectrum Mutants of c-KIT Kinase for Gastrointestinal Stromal Tumors. *J Med Chem.* 2019 Jul 11;62(13):6083-6101.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481