

1,6,8-trihydroxy-3-octyl-3,4-dihydro-1H-2-benzopyran-7-carboxylic acid

Chemical Properties


CAS No. :

Formula: C₁₈H₂₆O₆

Molecular Weight: 338.4

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,6,8-trihydroxy-3-octyl-3,4-dihydro-1H-2-benzopyran-7-carboxylic acid is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9551 mL	14.7754 mL	29.5508 mL
5 mM	0.591 mL	2.9551 mL	5.9102 mL
10 mM	0.2955 mL	1.4775 mL	2.9551 mL
50 mM	0.0591 mL	0.2955 mL	0.591 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481