


AI 3-25488

## Chemical Properties

CAS No. : 307-70-0  
 Formula: C<sub>11</sub>H<sub>4</sub>F<sub>20</sub>O  
 Molecular Weight: 532.12  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	AI 3-25488 is an agent of biochemical.
Targets(IC50)	Others

## Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8793 mL	9.3964 mL	18.7928 mL
5 mM	0.3759 mL	1.8793 mL	3.7586 mL
10 mM	0.1879 mL	0.9396 mL	1.8793 mL
50 mM	0.0376 mL	0.1879 mL	0.3759 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481