

**(±)-Abscisic acid****Chemical Properties**

CAS No. : 14375-45-2

Formula: C<sub>15</sub>H<sub>20</sub>O<sub>4</sub>


Molecular Weight: 264.32

Storage:

Keep away from direct sunlight, Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	(±)-Abscisic acid is an orally active phytohormone also found in animals. It regulates blood glucose in mammals.
Targets(IC50)	Others

**Solubility Information**

Solubility	DMSO: 240 mg/mL (907.99 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (7.57 mM), Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

**Preparing Stock Solutions**

	1mg	5mg	10mg
1 mM	3.7833 mL	18.9165 mL	37.8329 mL
5 mM	0.7567 mL	3.7833 mL	7.5666 mL
10 mM	0.3783 mL	1.8916 mL	3.7833 mL
50 mM	0.0757 mL	0.3783 mL	0.7567 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Wen F, et al. J Integr Plant Biol. 2015 Feb;57(2):213-28.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481