


UK-356202

Chemical Properties

CAS No. : 223671-94-1
 Formula: C₁₇H₁₃ClN₄O₂
 Molecular Weight: 340.76
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	UK-356,202 is a potent and selective urokinase-type plasminogen activator (K _i = 37 nM). Urokinase plasminogen activator (uPA, urokinase) is a trypsinlike serine protease and a therapeutic target for many cancer types, including breast, ovarian, and pancreatic cancer. uPA is a valuable oncology target.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9346 mL	14.6731 mL	29.3462 mL
5 mM	0.5869 mL	2.9346 mL	5.8692 mL
10 mM	0.2935 mL	1.4673 mL	2.9346 mL
50 mM	0.0587 mL	0.2935 mL	0.5869 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481