

## D-Leu-Ty hydrate

## Chemical Properties

CAS No. : 3303-29-5


Formula: C<sub>15</sub>H<sub>22</sub>N<sub>2</sub>O<sub>4</sub>

Molecular Weight: 294.35

Keep away from moisture, Store at low temperature

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	D-Leu-Ty hydrate ((S)-2-((R)-2-Amino-4-methylpentanamido)-3-(4-hydroxyphenyl) propanoic acid hydrate) is a tyrosine derivative commonly used in drug synthesis and biochemistry research.
Targets(IC50)	Amino Acids and Derivatives

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3973 mL	16.9866 mL	33.9732 mL
5 mM	0.6795 mL	3.3973 mL	6.7946 mL
10 mM	0.3397 mL	1.6987 mL	3.3973 mL
50 mM	0.0679 mL	0.3397 mL	0.6795 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

Luckose F, et al. Effects of amino acid derivatives on physical, mental, and physiological activities. Crit Rev Food Sci Nutr. 2015;55(13):1793-1144.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481