


GT-2227

Chemical Properties

CAS No. : 186096-20-8
 Formula: C₁₅H₂₄N₂
 Molecular Weight: 232.36
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	GT-2227 is a histamine H3 receptor antagonist.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3037 mL	21.5183 mL	43.0367 mL
5 mM	0.8607 mL	4.3037 mL	8.6073 mL
10 mM	0.4304 mL	2.1518 mL	4.3037 mL
50 mM	0.0861 mL	0.4304 mL	0.8607 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481