

## N-Mal-N-bis(PEG2-C2-Boc)

## Chemical Properties

CAS No. : 2100306-51-0


Formula: C<sub>29</sub>H<sub>48</sub>N<sub>2</sub>O<sub>11</sub>

Molecular Weight: 600.7

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	N-Mal-N-bis(PEG2-C2-Boc) is a polyethylene glycol (PEG)-based linker used in proteolysis targeting chimeras (PROTACs)[1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs, comprising two ligands joined by a linker—one targeting an E3 ubiquitin ligase and the other the target protein—leverage the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6647 mL	8.3236 mL	16.6472 mL
5 mM	0.3329 mL	1.6647 mL	3.3294 mL
10 mM	0.1665 mL	0.8324 mL	1.6647 mL
50 mM	0.0333 mL	0.1665 mL	0.3329 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

## Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481