

3-(4-Hydroxy-3-methoxyphenyl)propyl tetracosanoate

Chemical Properties


CAS No. :

Formula: C34H60O4

Molecular Weight: 532.85

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-(4-Hydroxy-3-methoxyphenyl)propyl tetracosanoate is a useful organic compound for research related to life sciences and the catalog number is T123820.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8767 mL	9.3835 mL	18.767 mL
5 mM	0.3753 mL	1.8767 mL	3.7534 mL
10 mM	0.1877 mL	0.9384 mL	1.8767 mL
50 mM	0.0375 mL	0.1877 mL	0.3753 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481