

Bis-propargyl-O-PEG17

Chemical Properties

CAS No. :


Formula: C₈H₁₀O₂

Molecular Weight: 138.16

Storage: Keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



17A

Biological Description

Description	Bis-propargyl-O-PEG17 is a PEG-based linker for PROTACs that connects two essential ligands, crucial for forming PROTAC molecules, facilitating selective protein degradation through the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are composed of two distinct ligands joined by a linker; one ligand targets an E3 ubiquitin ligase, and the other targets the desired protein. They utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	7.238 mL	36.1899 mL	72.3798 mL
5 mM	1.4476 mL	7.238 mL	14.476 mL
10 mM	0.7238 mL	3.619 mL	7.238 mL
50 mM	0.1448 mL	0.7238 mL	1.4476 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481