


Vitedoin A

Chemical Properties

CAS No. :	819861-40-0
Formula:	C ₂₀ H ₂₀ O ₆
Molecular Weight:	356.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Vitedoin A shows stronger antioxidative activity than alpha-tocopherol using the ferric thiocyanate method, it also shows higher radical-scavenging effect on the stable free radical, 1,1-diphenyl-2-picrylhydrazyl, than L-cysteine.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8061 mL	14.0304 mL	28.0607 mL
5 mM	0.5612 mL	2.8061 mL	5.6121 mL
10 mM	0.2806 mL	1.403 mL	2.8061 mL
50 mM	0.0561 mL	0.2806 mL	0.5612 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481