


Y 9000

Chemical Properties

CAS No. : 53944-40-4
 Formula: C₁₄H₉NO₃
 Molecular Weight: 239.23
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Y 9000 has anti-allergic activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.1801 mL	20.9004 mL	41.8008 mL
5 mM	0.836 mL	4.1801 mL	8.3602 mL
10 mM	0.418 mL	2.090 mL	4.1801 mL
50 mM	0.0836 mL	0.418 mL	0.836 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Maruyama Y, Terasawa M, Goto K, Oe T. [Anti-allergic activity of 7-acetyl-5-oxo-5H-[1] benzopyrano (2,3-b) pyridine (Y-9000) (author's transl)]. Nihon Yakurigaku Zasshi. 1978 Mar;74(2):179-91. Japanese. PubMed PMID: 77826.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481