

CP-467688

Chemical Properties


CAS No. : 909576-25-6

Formula: C₂₇H₂₇ClF₃N₃O₃

Molecular Weight: 533.97

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CP-467688 is a bio-active chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8728 mL	9.3638 mL	18.7276 mL
5 mM	0.3746 mL	1.8728 mL	3.7455 mL
10 mM	0.1873 mL	0.9364 mL	1.8728 mL
50 mM	0.0375 mL	0.1873 mL	0.3746 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481