

N-(Boc-PEG3)-N-bis(PEG3-acid)

Chemical Properties

CAS No. : 2055042-61-8


Formula: C₃₁H₆₀N₂O₁₅

Molecular Weight: 700.81

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	N-(Boc-PEG3)-N-bis(PEG3-acid) is a polyethylene glycol (PEG)-based linker compound used in the synthesis of proteolysis-targeting chimeras (PROTACs) [1].
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs consist of two ligands linked together: one binds to an E3 ubiquitin ligase and the other to the target protein, utilizing the intracellular ubiquitin-proteasome system for selective target protein degradation[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.4269 mL	7.1346 mL	14.2692 mL
5 mM	0.2854 mL	1.4269 mL	2.8538 mL
10 mM	0.1427 mL	0.7135 mL	1.4269 mL
50 mM	0.0285 mL	0.1427 mL	0.2854 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481