

Compound N040-0008

Chemical Properties


CAS No. :

Formula: C₄₆H₇₆O₁₇

Molecular Weight: 901.097

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N040-0008 is a useful organic compound for research related to life sciences and the catalog number is T131261.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1098 mL	5.5488 mL	11.0975 mL
5 mM	0.222 mL	1.1098 mL	2.2195 mL
10 mM	0.111 mL	0.5549 mL	1.1098 mL
50 mM	0.0222 mL	0.111 mL	0.222 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481