

Azide-PEG5-Boc

Chemical Properties

CAS No. : 1415800-41-7


Formula: C17H33N3O7

Molecular Weight: 391.46

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Azide-PEG5-Boc is a PEG-based linker for PROTACs that connects two essential ligands, facilitating the formation of PROTAC molecules. This linker enables selective protein degradation by leveraging the ubiquitin-proteasome system within cells.
Targets(IC50)	Others,PROTAC Linker
In vitro	PROTACs are chemical compounds featuring two distinct ligands linked by a connector; one ligand binds to an E3 ubiquitin ligase, while the other targets a specific protein. These compounds utilize the intracellular ubiquitin-proteasome system to selectively degrade target proteins[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5545 mL	12.7727 mL	25.5454 mL
5 mM	0.5109 mL	2.5545 mL	5.1091 mL
10 mM	0.2555 mL	1.2773 mL	2.5545 mL
50 mM	0.0511 mL	0.2555 mL	0.5109 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

An S, et al. Small-molecule PROTACs: An emerging and promising approach for the development of targeted therapy drugs. EBioMedicine. 2018 Oct;36:553-562

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481