


## Ganoderic acid AM1

## Chemical Properties

CAS No. : 149507-55-1  
 Formula: C<sub>30</sub>H<sub>42</sub>O<sub>7</sub>  
 Molecular Weight: 514.659  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
 Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Ganoderic acid AM1 treatment for 48 h inhibited the proliferation of HeLa human cervical carcinoma cells with IC(50) values of 19.8+/-0.7 microM.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.943 mL	9.7152 mL	19.4303 mL
5 mM	0.3886 mL	1.943 mL	3.8861 mL
10 mM	0.1943 mL	0.9715 mL	1.943 mL
50 mM	0.0389 mL	0.1943 mL	0.3886 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481