

Compound N050-0018

Chemical Properties


CAS No. :

Formula: C₂₃H₃₀O₆

Molecular Weight: 402.487

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound N050-0018 is a useful organic compound for research related to life sciences and the catalog number is T131250.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4845 mL	12.4227 mL	24.8453 mL
5 mM	0.4969 mL	2.4845 mL	4.9691 mL
10 mM	0.2485 mL	1.2423 mL	2.4845 mL
50 mM	0.0497 mL	0.2485 mL	0.4969 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481