


2''-O-Acetylsprengerinin C

Chemical Properties

CAS No. :	1220707-33-4
Formula:	C ₄₆ H ₇₂ O ₁₇
Molecular Weight:	897.065
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	2''-O-Acetylsprengerinin C is a natural product for research related to life sciences. The catalog number is TN2782 and the CAS number is 1220707-33-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1147 mL	5.5737 mL	11.1474 mL
5 mM	0.2229 mL	1.1147 mL	2.2295 mL
10 mM	0.1115 mL	0.5574 mL	1.1147 mL
50 mM	0.0223 mL	0.1115 mL	0.2229 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481