

## CXCR2 antagonist 5

## Chemical Properties

CAS No. :

Formula: C15H14F2N4O2S

Molecular Weight: 352.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	CXCR2 antagonist 5 (compound 25) is a potent CXCR2 antagonist with strong binding affinity (IC50 = 0.013 $\mu$ M) and effective calcium mobilization (IC50 = 0.1 $\mu$ M) [1].
Targets(IC50)	Others,CXCR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.838 mL	14.190 mL	28.3801 mL
5 mM	0.5676 mL	2.838 mL	5.676 mL
10 mM	0.2838 mL	1.419 mL	2.838 mL
50 mM	0.0568 mL	0.2838 mL	0.5676 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481