


Hexbutinol methiodide

Chemical Properties

CAS No. :	127471-24-3
Formula:	C ₂₂ H ₃₄ INO ₂
Molecular Weight:	471.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Hexbutinol methiodide is an inhibitor of muscarinic receptors.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1213 mL	10.6063 mL	21.2125 mL
5 mM	0.4243 mL	2.1213 mL	4.2425 mL
10 mM	0.2121 mL	1.0606 mL	2.1213 mL
50 mM	0.0424 mL	0.2121 mL	0.4243 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481