

## New compound 3

## Chemical Properties

CAS No. :

Formula: C<sub>27</sub>H<sub>30</sub>O<sub>15</sub>

Molecular Weight: 594.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description

New compound 3 can be used for related research in the field of life sciences, and its product number is TN6410.

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6821 mL	8.4104 mL	16.8209 mL
5 mM	0.3364 mL	1.6821 mL	3.3642 mL
10 mM	0.1682 mL	0.841 mL	1.6821 mL
50 mM	0.0336 mL	0.1682 mL	0.3364 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481