

m-PEG8-Ms

Chemical Properties

CAS No. : 1059588-19-0


Formula: C18H38O11S

Molecular Weight: 462.55

Store at low temperature

Storage: Pure form: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	m-PEG8-Ms is a PROTAC linker belonging to the PEG class and can be used to synthesize PROTAC compounds.
Targets(IC50)	ADC Linker, PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1619 mL	10.8096 mL	21.6193 mL
5 mM	0.4324 mL	2.1619 mL	4.3239 mL
10 mM	0.2162 mL	1.081 mL	2.1619 mL
50 mM	0.0432 mL	0.2162 mL	0.4324 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sanxing Sun, et al. Triazolotriazine derivatives as a2a receptor antagonists. WO2020002969A1.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481