


3-Deoxy-1,2-O-isopropylidene-5-O-(p-toluoyl)- α -D-glycero-pent-3-enofuranose

Chemical Properties

CAS No. :	75096-63-8
Formula:	C ₁₆ H ₁₈ O ₅
Molecular Weight:	290.31
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-Deoxy-1,2-O-isopropylidene-5-O-(p-toluoyl)- α -D-glycero-pent-3-enofuranose is a valuable organic compound for life sciences research. Its catalog number is TNU1040 and its CAS number is 75096-63-8.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4446 mL	17.223 mL	34.4459 mL
5 mM	0.6889 mL	3.4446 mL	6.8892 mL
10 mM	0.3445 mL	1.7223 mL	3.4446 mL
50 mM	0.0689 mL	0.3445 mL	0.6889 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481