

Aniline, p-(5-phenoxypropyloxy)-

Chemical Properties


CAS No. : 101582-06-3

Formula: C₁₇H₂₁NO₂

Molecular Weight: 271.35

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Aniline, p-(5-phenoxypropyloxy)- is a bioactive chemical.
Targets(IC ₅₀)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6853 mL	18.4264 mL	36.8528 mL
5 mM	0.7371 mL	3.6853 mL	7.3706 mL
10 mM	0.3685 mL	1.8426 mL	3.6853 mL
50 mM	0.0737 mL	0.3685 mL	0.7371 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481