


## Dammaradienyl acetate

## Chemical Properties

CAS No. :	52914-31-5
Formula:	C <sub>32</sub> H <sub>52</sub> O <sub>2</sub>
Molecular Weight:	468.766
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Dammaradienyl acetate is a natural product of Inula, Asteraceae. The catalog number is TN3756 and the CAS number is 52914-31-5. Dammaradienyl acetate can be used as a reference standard.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1332 mL	10.6662 mL	21.3324 mL
5 mM	0.4266 mL	2.1332 mL	4.2665 mL
10 mM	0.2133 mL	1.0666 mL	2.1332 mL
50 mM	0.0427 mL	0.2133 mL	0.4266 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481