

Lucidenic acid F

Chemical Properties


CAS No. : 98665-18-0

Formula: C₂₇H₃₆O₆

Molecular Weight: 456.57

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Lucidenic acid F as a modulator of JNK and p38, it shows potent inhibitory effects on EBV-EA induction.
Targets(IC50)	JNK,p38 MAPK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1902 mL	10.9512 mL	21.9024 mL
5 mM	0.438 mL	2.1902 mL	4.3805 mL
10 mM	0.219 mL	1.0951 mL	2.1902 mL
50 mM	0.0438 mL	0.219 mL	0.438 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481