

AA-dUTP sodium salt

Chemical Properties

CAS No. : 936327-10-5


Formula: C₁₂H₂₀N₃O₁₄P₃.xNa

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	AA-dUTP [sodium salt] is a fluorescent dye used for staining cDNA.
Targets(IC50)	Others,DNA/RNA Synthesis

Reference

Cox WG, et al. Fluorescent DNA hybridization probe preparation using amine modification and reactive dye coupling. *Biotechniques*. 2004 Jan;36(1):114-22.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481