


Ciproquazone

Chemical Properties

CAS No. :	33453-23-5
Formula:	C ₁₉ H ₁₈ N ₂ O ₂
Molecular Weight:	306.36
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Ciproquazone has analgesic activity and is used to treat spondylitis.
Targets(IC50)	Others, Prostaglandin Receptor

Solubility Information

Solubility	DMSO: 60 mg/mL (195.85 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2641 mL	16.3207 mL	32.6413 mL
5 mM	0.6528 mL	3.2641 mL	6.5283 mL
10 mM	0.3264 mL	1.6321 mL	3.2641 mL
50 mM	0.0653 mL	0.3264 mL	0.6528 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Vega S, et al. Thiophene isosteres: synthesis and biological evaluation of 3-substituted derivatives of 4-phenyl-2-thioxo-benzo [4, 5] thieno [2, 3-d] pyrimidine. European journal of medicinal chemistry, 1991, 26(3): 323-329.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481